HydroCision®

The Leader in Hydrosurgery®

HydroCision SpineJet® Discectomy Devices

Symbol Description

	Catalogue Number		Do Not Re-Sterilize
	Lot Number		Operating Range
	Do Not Use If Package Is Damaged		Expiration Date
	Single Use		Sterilize Ethylene Oxide
	Manufacturer		Consult for Use
	Authorized Representative in the European Community		U.S Federal law restricts this device to sales by or on the order of a physician.
	Connection of the waste hose, or any container connected to it, to a vacuum source is not required and will create continuous suction at the device tip and affect device performance		

English

Instructions for Use

Please read all information carefully. Also refer to the Product Manual supplied with the HydroCision SpineJet Console. The components of the HydroCision SpineJet Discectomy Devices are designed to be used in combination with the HydroCision SpineJet Console and function as a single unit. Failure to properly follow instructions may lead to mechanical, electrical, or thermal injury and result in improper functioning of the device.

HydroCision, Inc.
267 Boston Road, Suite 28
North Billerica, MA 01862
USA
Customer Care:
(978) 474-9300

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 (0) 70 345 8570

Indications

The HydroCision SpineJet Discectomy Devices are indicated for the cutting and removal of soft tissue and the ablation and removal of hard tissue or bone is required. Specific functions include cutting, ablation and shaping of soft tissue, and decorticating and smoothing of bone, cartilage and other bone related tissue in open and minimally invasive disc decompression and spinal fusion surgeries.

Disposable Assembly Description

The Discectomy Devices consists of a disposable hand piece that is connected to the console through a disposable hose and a pump cartridge.

The Pump Cartridge mounts into a nest assembly on the front of the HydroCision SpineJet Console; this connection provides the power to the Discectomy Devices. A fluid supply hose connects via a standard bag spike to the sterile irrigant bag. A waste hose from the Handpiece connects via a standard connector to a waste container.

The SpineJet Percutaneous Discectomy and MicroDiscectomy include respective disposable and reusable stainless steel guide wire, spinal needle, introducer, and dilator. The Access Set is used to provide a surgical pathway for the SpineJet Discectomy procedures.

Warnings

- All HydroCision SpineJet Discectomy Devices should not be activated in close proximity to, or come into contact with, the spinal cord, nerve roots, or major blood vessels to avoid the possibility of injury or damage from the fluidjet cutting action, locally generated vacuum, or device edges or tips.
- Do not touch the SpineJet hand piece tip while activated.
- Do not insert or withdraw the SpineJet hand piece while activated.
- Inadvertent activation or movement of the SpineJet hand piece or Access Set outside the field of vision or without adequate assurance of device placement via fluoroscopy or an alternate imaging technology may result in patient injury.
- Do not reuse any system component or accessory labeled as SINGLE USE. Reuse may lead to infection or injury to the patient.
- Attempts to bend the SpineJet hand piece or Access Sets may render the tool unusable or unsafe.
- SpineJet Devices utilizing an Access Introducer Set should be removed through the access introducer prior to removing the access introducer from the surgical site.
- Depth markings are provided on the instrument shaft to assist in device placement: the spaced band markings serve as visual indicators to avoid unintended puncture of the anterior annulus

and to alert the user to the proximity of the fluidjet to the annulotomy.

- No modification of this equipment is allowed

Precautions

- A thorough understanding of the principles and techniques involved in spinal procedures is essential to avoid injury to the patient and medical personnel, and damage to the device or other medical instruments.
- Read all instructions carefully. Failure to properly follow instructions may lead to electrical, mechanical, or thermal injury and cause improper functioning of the device.
- Please refer to the HydroCision SpineJet Console Product Manual for step-by-step instructions regarding the assembly and initial system check of the HydroCision SpineJet Console.
- SpineJet Disposable Handpiece Assemblies and Percutaneous Discectomy Access Sets are supplied as sterile, single use disposable devices. If the package is opened or damaged the sterility of the hand piece will be compromised.
- Ensure complete connection of the Handpiece assembly to the Console, fluid supply and waste container prior to use.
- Avoid kinking of the high-pressure hose, which can cause leaks.
- The SpineJet hand piece or Access Set should be inserted, manipulated, and withdrawn carefully from the operative site to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Use of the higher settings on the device console will lead to more aggressive tissue removal. Use caution near sensitive tissues, such as neurovascular bundles.
- It is recommended that the SpineJet Disposable Assembly be used on console setting #10 for nucleus removal.
- This device can cut soft tissue.
- A transdural approach should not be used under any circumstances.
- Do not apply excessive force in any direction during the procedure to avoid patient injury.
- Particular precaution should be made to avoid unintended puncture of the annulus.
- MicroResector and EndoResector hand pieces are provided with an o-ring depth indicator used to set the desired working length beyond the distal end of a scope. Prior to surgery, set the desired working length by carefully advancing the Resector through the working channel of the scope and sliding the o-ring into position.
- The SpineJet Access Sets are not intended for injection.

Instructions for Use

Also see the HydroCision Console Product Manual.

1. Circulating Nurse: Ensure that the sealed package is undamaged. Using sterile technique, carefully open the product outer package and present contents to sterile field personnel.

2. Scrub Nurse or Surgeon: Open inner pouch or lid.
3. Scrub Nurse or Surgeon: Clip the high-pressure hose/waste hose to the sterile drape.
4. Scrub Nurse or Surgeon: Pass the pump cartridge, high-pressure hose/waste hose, and coiled supply hose with bag spike to the Circulating Nurse.
5. Circulating Nurse: Open the Console nest by turning the knob to the right. Insert the pump cartridge into the nest by pushing the cartridge in until it is fully seated, and close the nest by turning the knob to the left. Attach the waste hose connector to a waste collection container – connection of the waste collection container to a vacuum source is not required and will create continuous suction at the device tip and affect device performance. The waste container must be vented to the atmosphere if not connected to a vacuum source. Remove the sterile cover from the bag spike and insert into an irrigant supply bag. A 1-liter bag is recommended for discectomies. Ensure that there are no kinks or external obstructions in the supply, high pressure, or waste hoses.
6. Circulating Nurse: Ensure that the foot switch is plugged into the connector on the front of the console and the power cord is plugged into the back of the console. Plug the power cord into a proper electrical outlet. Turn on the main power switch and check that the green ready light is illuminated.
7. Nursing Staff: Once the system has been primed with saline, do not allow the saline bag to empty; an empty bag will allow air into the system and reduce the system's efficiency. Therefore, always change to a new saline bag before the bag in use empties. Take care when switching bags to prevent air from entering system by closing off the supply hose with the pinch clamp.
8. See Percutaneous Access Set if using SpineJet Percutaneous HydroDiscectomy Set. See MicroResector Access Set if using MicroResector Set.
9. Prepare the patient pre-operatively according to standard procedures.
10. It is recommended that the SpineJet be used on console setting #10 for nucleus removal.
11. After completing the procedure, disconnect the SpineJet Disposable Assembly from the Power Console by turning the nest lever to the right and removing the pump cartridge by pulling it straight out.

Percutaneous Access Set:

1. Carefully remove Access Set tray from the inner pouch while being careful of all sharp edges. Check to ensure the dilator and spinal needle tips are undamaged.
2. Fully insert the polymer dilator into the selected introducer while locking the hub in place with a ¼ turn to the right. Assure

- proper locking by inspecting the distal tip for transition alignment.
3. The SpineJet PercResector and all Access Set components should be used under local anesthesia for percutaneously accessed procedures to allow for patient monitoring for signs of nerve root irritation.
 4. The included spinal needle should be advanced from a standard extra-pedicular posterior-lateral approach inferior to the neural foramen into the center of the disc using a standard discography approach under fluoroscopic guidance. A transdural approach should not be used under any circumstances.
 5. Once this approach has been completed for the target disc level, remove the stylet from the spinal needle. Care must be taken to ensure that the spinal needle remains in place.
 6. The Percutaneous Access Set guide wire should be advanced through the spinal needle into the nucleus while ensuring the spinal needle remains in place.
 7. Confirm proper guide wire position with A/P and lateral fluoroscopy views. Dilute contrast media may help identify the nucleus boundary. Once proper placement is confirmed, remove the spinal needle by sliding it up and over the guide wire. Care must be taken to ensure that the guide wire remains in place.
 8. While holding the guide wire in place slide the dilator/introducer as a unit over the guide wire and insert the dilator/introducer unit through the annulus into the nucleus using a twisting motion. Care should be taken to ensure that the dilator remains fully inserted into the introducer until the nucleus is reached.
 9. Confirm position with fluoroscopy and remove the guide wire and dilator. Care must be taken to ensure that the introducer remains in place.
 10. Insert and run the SpineJet Percutaneous Discectomy Device. Do not apply excessive force in any direction during the procedure to avoid patient injury.
 11. When using the SpineJet Discectomy Devices stop the procedure if patient complains of sudden onset of pain.
 12. Remove the SpineJet Discectomy Device's instrument shaft from the patient prior to removing the introducer. Follow standard surgical procedure for post-operative cleaning and closure of the percutaneous site.
 13. If the device tip becomes blocked with foreign matter, it will typically be noticed by a reduction in device efficiency or the presence of spray from the tip; b) Stop the jet flow by releasing the foot switch; c) Remove the hand piece from the surgical site using care not to come into contact with vital structures; d) Remove the obstruction with forceps taking care not to touch the opening in the high-pressure jet. Once removed, depress the foot switch and check that there is a single coherent jet flow. If the obstruction is not completely removed, repeat procedure.

MicroResector Access Set: (Sterilize before Use)

1. Also refer to the Instructions for Use supplied with the HydroCision SpineJet MicroResector Disposable Assembly.
2. Prior to the surgical procedure, clean and sterilize the SpineJet MicroDiscectomy Access Handpiece Set through a facility-approved process for stainless steel surgical instruments. Particular care must be taken to assure that the lumens of both the handpiece and dilator are free of debris prior to sterilization.
3. Inspect for sharpness, wear, or damage prior to every use.
4. Ensure that the needle will pass freely through the dilator and that the dilator will pass freely through the handpiece prior to use.
5. The Access Set needle or solid one step dilator should be advanced into the annulus adjacent to either the traversing or exiting nerve root. The dilator threaded knob should be attached for easy insertion.
6. Confirm proper instrument position (needle or dilator) before proceeding with handpiece placement with fluoroscopy views. If using the needle slide the dilator over the needle (if used) and insert the dilator through the annulus into the nucleus.
7. Remove the dilator threaded knob and slide the introducer handpiece over the dilator through the annulus and into the nucleus. Remove the needle if used. Note that on some handpieces there are two markings on the handpiece. The distal mark is not a measurement and is designed to signal that the end of the handpiece is ready to come out of the annulotomy. The 1cm depth marking on the handpiece can be used to judge initial handpiece depth; however fluoroscopy should be used to confirm final position. Other handpieces have a tactile indicator when the handpiece is inserted 7.5 mm into the annulus.
8. Insert and run the SpineJet MicroResector. Use a gentle pistoning motion while slightly rotating and angulating the instrument shaft to consume nucleus. Adequate nucleus consumption is typically achieved in 2-3 minutes run time. Do not apply excessive force in any direction during the procedure to avoid patient injury.
9. Remove the SpineJet MicroResector instrument shaft from the patient prior to removing the handpiece.

After Surgery

The SpineJet hand piece assembly, Saline Bag, and waste receptacle collection container may be discarded using biohazard disposal procedures according to local regulations.

How Supplied

The SpineJet handpiece assembly is provided sterile. The contents are sterile unless the package is opened or damaged. Do not resterilize. Do not use if package is opened or damaged.

Storage

The SpineJet® Discectomy Devices should be stored in a manner to prevent damage to the instruments prior to use.

Temperature Range: -40°F (-40°C) to 125°F (52°C)

Humidity Range: 0% to 100%, noncondensing


Atmospheric Pressure: 500 to 1500 millibar

For further information, consult the HydroCision Console Product Manual.

HydroCision, SpineJet, and The Leader in Hydrosurgery are registered trademarks of HydroCision, Inc. ©2004-2006. All rights reserved.

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