HydroCision®

The Leader in Hydrosurgery®

HydroCision SpineJet® Disposable Assembly with Quick Connect

Symbol Description


Catalogue Number	Do Not Re-Sterilize	Lot Number
Do Not Use If Package Is Damaged	Operating Range	Expiration Date
Manufacturer	Sterilized Using Ethylene Oxide	Single Use
Consult Instructions for Use	Authorized Representative in the European Community	
U.S. Federal law restricts this device to sales by or on the order of a physician.		
Connection of the waste hose, or any container connected to it, to a vacuum source is not required and will create continuous suction at the device tip and affect device performance		

Instructions for Use

Please read all information carefully. Also refer to the Product Manual supplied with the HydroCision SpineJet Console. The components of the HydroCision SpineJet Disposable Assembly are designed to be used in combination with the HydroCision SpineJet Console and function as a single unit. Failure to properly follow instructions may lead to mechanical, electrical, or thermal injury and result in improper functioning of the device.

HydroCision, Inc.
267 Boston Road, Suite 28
North Billerica, MA 01862 USA
Customer Care: (978) 474-9300

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 (0) 70 345 8570

Indications

The HydroCision SpineJet System is indicated for surgical procedures where the cutting and removal of soft tissue and the ablation and removal of hard tissue or bone is required. Specific functions include cutting, ablation and shaping of soft tissue, and decorticating and smoothing of bone, cartilage and other bone related tissue in open and minimally invasive disc decompression and spinal fusion surgeries.

Disposable Assembly Description

The SpineJet Disposable Assembly consists of a disposable handpiece that is connected to the console through a disposable Quick Connect, a Pump Cartridge, and hoses. The Pump Cartridge mounts into a nest assembly on the front of the HydroCision SpineJet Console; this connection provides the power to the Disposable Assembly. A fluid supply hose connects via a standard bag spike to the sterile irrigant bag. A waste hose from the Quick Connect connects via a standard connector to a waste container.

Warnings

- The HydroCision SpineJet Disposable Assembly should not be activated in close proximity to, or come into contact with, the spinal cord, nerve roots, or major blood vessels to avoid the possibility of injury or damage from the fluidjet cutting action, locally generated vacuum, or device edges or tip.
- Do not touch the SpineJet handpiece tip while activated.
- Do not insert or withdraw the SpineJet handpiece while activated.
- Inadvertent activation or movement of the SpineJet handpiece outside the field of vision or without adequate assurance of device placement via fluoroscopy or an alternate imaging technology may result in patient injury.
- Do not reuse any system component or accessory labeled as SINGLE USE. Reuse may lead to infection or injury to the patient.
- Attempts to bend the SpineJet handpiece may render the tool unusable or unsafe.
- SpineJet MicroResector and SpineJet MicroResector MD instruments should be removed through the access cannula prior to removing the access cannula from the surgical site.
- Depth markings are provided on the instrument shaft to assist in device placement: the spaced band

markings serve as visual indicators to avoid unintended puncture of the anterior annulus and to alert the user to the proximity of the fluidjet to the annulotomy.

- No modification of this equipment is allowed

Precautions

- A thorough understanding of the principles and techniques involved in spinal surgeries is essential to avoid injury to the patient and medical personnel, and damage to the device or other medical instruments.
- Read all instructions carefully. Failure to properly follow instructions may lead to electrical, mechanical, or thermal injury and cause improper functioning of the device.
- Please refer to the HydroCision SpineJet Console Product Manual for step-by-step instructions regarding the assembly and initial system check of the HydroCision SpineJet Console.
- SpineJet Disposable Assemblies and Quick Connects are supplied sterile. If the package is opened or damaged the sterility of the handpiece will be compromised.
- Ensure complete connection of SpineJet handpiece to the Quick Connect prior to use.
- Ensure complete connection of the Quick Connect assembly to the Console, fluid supply and waste container prior to use.
- Avoid kinking of the high-pressure hose, which can cause leaks.
- The SpineJet handpiece should be inserted, manipulated, and withdrawn carefully from the operative site to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Use of the higher settings on the device console will lead to more aggressive tissue removal. Use caution near sensitive tissues, such as neurovascular bundles.
- It is recommended that the SpineJet Disposable Assembly be used on console setting #10 for nucleus removal.
- This device can cut soft tissue.
- A transdural approach should not be used under any circumstances.
- Do not apply excessive force in any direction during the procedure to avoid patient injury.
- Particular precaution should be made to avoid unintended puncture of the annulus.

- To prevent clogging of the device tip, avoid applying excessive force with the cutting and scraping edges of the XL devices that would release cartilage fragments larger than the device's evacuation tube opening.
- The EndoResector handpiece is provided with an o-ring depth indicator used to set the desired working length beyond the distal end of an endoscope. Prior to surgery, set the desired working length by carefully advancing the EndoResector through the working channel of the endoscope and sliding the o-ring into position.

Instructions for Use

Also see the HydroCision Console Product Manual.

1. Circulating Nurse: Ensure that the sealed package is undamaged. Using sterile technique, carefully open the product outer package and present contents to sterile field personnel.
2. Scrub Nurse or Surgeon: Open inner pouch or lid. If present, remove the plastic guard from the distal tip of the handpiece.
3. Attach the SpineJet Disposable Assembly to the Quick Connect by aligning the alignment features on the handpiece handle with corresponding features on the Quick Connect and sliding them into place until the locking tabs click indicating the correct attachment. Do not force the connection. Note that when energized, the Quick Connect will slide away from the handpiece, seating it into the locking position. This creates a small gap between the Quick Connect and handpiece.
4. Scrub Nurse or Surgeon: Clip the high-pressure hose/waste hose to the sterile drape.
5. Scrub Nurse or Surgeon: Pass the pump cartridge, high-pressure hose/waste hose, and coiled supply hose with bag spike to the Circulating Nurse.
6. Circulating Nurse: Open the Console nest by turning the knob to the right. Insert the pump cartridge into the nest by pushing the cartridge in until it is fully seated, and close the nest by turning the knob to the left. Attach the waste hose connector to a waste collection container – connection of the waste collection container to a vacuum source is not required and will create continuous suction at the device tip and affect device performance. The waste container must be vented to the atmosphere if not connected to a vacuum source. Remove the sterile cover from the bag spike and insert into an irrigant supply bag. A 3 liter bag is recommended for fusion

procedures and a 1-liter bag is recommended for discectomies. Ensure that there are no kinks or external obstructions in the supply, high pressure, or waste hoses.

7. Circulating Nurse: Ensure that the foot switch is plugged into the connector on the front of the console and the power cord is plugged into the back of the console. Plug the power cord into a proper electrical outlet. Turn on the main power switch and check that the green ready light is illuminated.
 8. Nursing Staff: Once the system has been primed with saline, do not allow the saline bag to empty; an empty bag will allow air into the system and reduce the system's efficiency. Therefore, always change to a new saline bag before the bag in use empties. Take care when switching bags to prevent air from entering system by closing off the supply hose with the pinch clamp.
 9. a) If the device tip becomes blocked with foreign matter, it will typically be noticed by a reduction in device efficiency or the presence of spray from the tip; b) Stop the jet flow by releasing the foot switch; c) Remove the handpiece from the surgical site using care not to come into contact with vital structures; d) Remove the obstruction with forceps taking care not to touch the opening in the high-pressure jet. Once removed, depress the foot switch and check that there is a single coherent jet flow. If the obstruction is not completely removed, repeat procedure.
- To Change Tools:*
- a. Unseat the locking tabs by pushing the Quick Connect and handpiece completely together.
 - b. Depress the two locking tabs at opposite sides of the Quick Connect and gently remove the tool from the Quick Connect while keeping the tabs fully depressed.
 - c. Insert the next desired tool as above.
10. After completing the procedure, disconnect the SpineJet Disposable Assembly from the Power Console by turning the nest lever to the right and removing the pump cartridge by pulling it straight out. If desired, the SpineJet Disposable Assembly can be disconnected from the Quick Connect by depressing the two locking tabs at opposite sides on the proximal end of the tool and pulling the tool from the Quick Connect while keeping the tabs depressed.

After Surgery

The SpineJet handpiece assembly, Quick Connect, Saline Bag, and waste receptacle collection container may be discarded using biohazard disposal procedures according to local regulations.

How Supplied

The SpineJet handpiece assembly and Quick Connect are provided sterile.

The contents are sterile unless the package is opened or damaged. Do not resterilize.

Do not use if package is opened or damaged.

Storage


STERILE	EO
---------	----

Do not store product above 125°F (52°C) or below -40°F (-40°C). Avoid storage near moisture and direct heat.

For further information, consult the HydroCision Console Product Manual.

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