## HydroCision® SpineJet® Disposable Handpiece.

### Symbol Description

	Catalogue Number		Do Not Re-Sterilize
	Lot Number		Operating Range
	Do Not Use If Package Is Damaged		Expiration Date
	STERILIZED USING CONSULT		Single Use
	Manufacturer Instructions for Use		Sterilized Using Consult Ethylene Oxide

Authorized Representative in the European Community

U.S Federal law restricts this device to sales by or on the order of a physician.

Connection of the waste hose, or any container connected to it, to a vacuum source is not required and will create continuous suction at the device tip and affect device performance


### Instructions for Use

Please read all information carefully. Also refer to the Product Manual supplied with the HydroCision SpineJet console. The components of the HydroCision SpineJet Disposable handpiece is designed to be used in combination with the HydroCision SpineJet console and function as a single unit. Failure to properly follow instructions may lead to mechanical, electrical, or thermal injury and result in improper functioning of the device.

HydroCision, Inc.  
267 Boston Road, Suite 28  
North Billerica, MA 01862 USA  
Customer Care: (978) 474-9300

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
+31 (0) 70 345 8570

### Indications

The HydroCision SpineJet System is indicated for surgical procedures where the cutting and removal of soft tissue and the ablation and removal of hard tissue or bone is required. Specific functions include cutting, ablation and shaping of soft tissue, and decorticating and smoothing of bone, cartilage and other bone related tissue in open and minimally invasive in spinal surgeries.

### Disposable Device Description

The SpineJet Hydrosurgery System handpiece consists of a disposable handpiece that is connected to the console with a high pressure hose and pump cartridge. The Pump Cartridge mounts into a nest assembly on the front of the HydroCision Console; this connection provides the power to the handpiece. A fluid supply hose connects via a standard bag spike to the sterile irrigant bag. A waste hose connects via a standard connector to a waste container.

### Warning

- The SpineJet should not be activated in close proximity to, or come into contact with, the spinal cord, nerve roots, or major blood vessels to avoid the possibility of injury or damage from the fluidjet cutting action, locally generated vacuum, or device edges or tip.
- Do not touch the SpineJet tip while activated.
- Do not insert or withdraw the SpineJet while activated.
- Inadvertent activation or movement of the SpineJet outside the field of vision or without adequate assurance of device placement via fluoroscopy or an alternate imaging technology may result in patient injury. **Depth markings** are provided on the instrument shaft to assist in device placement: the distal marking is located 1cm from the evacuation tube edge and warns the user that the jet cutting action is about to be removed from the disc space. The proximal markings allow the user to insert the device 0.5 cm or 1.0 cm past the end of the cannula within the disc space. These band markings serve as visual indicators to avoid unintended puncture of the anterior annulus.
- Do not reuse any system component or accessory labeled as SINGLE USE. Reuse may lead to infection or injury to the patient.
- Attempts to bend the SpineJet may render the tool unusable or unsafe.
- The SpineJet instrument should be removed through the access cannula prior to removing the access cannula from the surgical site.
- No modification of this equipment is allowed

### Precautions

- A thorough understanding of the principles and techniques involved in spinal surgeries is essential to avoid injury to the patient and medical personnel, and damage to the device or other medical instruments.

- Read all instructions carefully. Failure to properly follow instructions may lead to electrical, mechanical, or thermal injury and cause improper functioning of the device.
- Please refer to the HydroCision Console Product Manual for step-by-step instructions regarding the assembly and initial system check of the HydroCision Console.
- SpineJet handpiece is supplied sterile. If the package is opened or damaged the sterility of the handpiece will be compromised.
- Ensure complete connection of SpineJet handpiece to the console, fluid supply and waste container prior to use.
- Avoid kinking of the high-pressure hose, which can cause leaks.
- The SpineJet should be inserted, manipulated, and withdrawn carefully from the operative site to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- A transdural approach should not be used under any circumstances.
- Do not apply excessive force in any direction during the procedure to avoid patient injury.
- Use of the higher settings on the device console will lead to more aggressive tissue removal.
- Use caution near sensitive tissues, such as neurovascular bundles.
- This device can cut soft tissue.
- Particular precaution should be made to avoid unintended puncture of the annulus.

### Instructions for Use

Also see the HydroCision Console Product Manual.

1. Circulating Nurse: Ensure that the sealed package is undamaged. Using sterile technique, carefully open the product packaging and present contents to sterile field personnel.
2. Scrub Nurse or Surgeon: Open inner pouch.
3. Scrub Nurse or Surgeon: Remove the disposable handpiece assembly from the inner sleeve packaging (see label for designation) and clip the high-pressure hose/waste hose to the sterile drape leaving enough length of hose to perform the procedure.
4. Scrub Nurse or Surgeon: Pass the inner sleeve packaging containing the console components to the Circulating Nurse.
5. Circulating Nurse: Open the Console nest by turning the knob to the right. Insert the pump cartridge into the nest by pushing the cartridge in until it is fully seated, and close the nest by turning the knob to the left. Attach the waste hose connector to a waste collection container – connection of the waste collection container to a vacuum source is not required and will create continuous suction at the device tip and affect device performance. The waste container must be vented to

the atmosphere if not connected to a vacuum source. Remove the sterile cover from the bag spike and insert into an irrigant supply bag. A 1-liter bag is recommended. Ensure that there are no kinks or external obstructions in the supply, high pressure, or waste hoses.

6. Circulating Nurse: Ensure that the foot switch is plugged into the connector on the front of the console and the power cord is plugged into the back of the console. Plug the power cord into a proper electrical outlet. Turn on the main power switch and check that the green ready light is illuminated.
7. Nursing Staff: Once the system has been primed with saline, do not allow the saline bag to empty; an empty bag will allow air into the system and reduce the system's efficiency. Therefore, always change to a new saline bag before the bag in use empties. Take care when switching bags to prevent air from entering system by closing off the supply hose with the pinch clamp.
8. Prepare the patient pre-operatively according to standard procedures. After making a standard laminotomy and identifying the nerve root, the SpineJet® Access Cannula is inserted into the annulus adjacent to the traversing nerve root. (refer to Access Cannula instructions for use p/n 1000-2348.).
9. The o-ring depth indicator is used to set the desired working length of the SpineJet from the end of the Access Cannula. The o-ring depth indicator is initially set and comes out of the sterile packaging set at a depth of 0.5 cm. If desired the surgeon may slide the o-ring to the most proximal shaft marking indication a depth of 1.0 cm. It is not recommended that the device be inserted any deeper than 1cm. Carefully advance the SpineJet instrument shaft into the introducer cannula and confirm proper position with fluoroscopy.
10. It is recommended that the SpineJet be used on console setting #10 for nucleus removal.
11. Activate the footswitch to begin the discectomy procedure. Use a gentle pistoning motion while slightly rotating and angulating the instrument shaft to consume nucleus. Adequate nucleus consumption is typically achieved in 2-3 minutes run time. Do not apply excessive force in any direction during the procedure to avoid patient injury.
12. If the device tip becomes blocked with foreign matter, it will typically be noticed by a reduction in device efficiency or the presence of spray from the tip. Stop the jet flow by releasing the foot switch and remove the obstruction with forceps taking care not to touch the opening in the high-pressure jet. Once removed, depress the foot switch and check that there is a single coherent jet flow. If the obstruction is not completely removed, repeat procedure.
13. Remove the SpineJet instrument shaft from the patient prior to removing the introducer cannula.
14. After completing the procedure, disconnect the SpineJet handpiece from the Power Console by turning the nest lever to the right and removing the pump cartridge by pulling it straight out.

#### After Surgery

The SpineJet handpiece, saline bag, and waste receptacle collection container may be discarded using standard biohazard disposal procedures according to local regulations.

#### How Supplied

The SpineJet handpiece is provided sterile.

The contents are sterile unless the package is opened or damaged.

Do not resterilize.


Do not use if package is opened or damaged.

#### Storage

Do not store product above 125°F (52°C) or below -40°F (-40°C).

Avoid storage near moisture and direct heat.

#### Additional Symbol Description

**Console Components:** Items within this container are intended to interface with the SpineJet Console, waste container, and sterile saline supply

**Patient Components:** The surgical handpiece is located within this container

For further information, consult the HydroCision® Console Product Manual.

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